

小鼠抗 GTF2B 单克隆抗体

中文名称: 小鼠抗 GTF2B 单克隆抗体

英文名称: Anti-GTF2B mouse monoclonal antibody

别名: general transcription factor IIB; TF2B; TFIIB

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Mouse

抗原: GTF2B

反应种属: Human, Mouse, Monkey

标记物: Unconjugate

克隆类型: Mouse Monoclonal

技术规格

Background:	This gene encodes the general transcription factor IIB, one of the ubiquitous factors required for transcription initiation by RNA polymerase II. The protein localizes to the nucleus where it forms a complex (the DAB complex) with transcription factors IID and IIA. Transcription factor IIB serves as a bridge between IID, the factor which initially recognizes the promoter sequence, and RNA polymerase II.
Applications:	WB
Name of antibody:	GTF2B
Immunogen:	Fusion protein of human GTF2B
Full name:	general transcription factor IIB

Synonyms:	TF2B; TFIIB
SwissProt:	Q00403
WB Predicted band size:	35 kDa
WB Positive control:	3T3, COS7, MCF7, Jurkat and Hela cell lysates
WB Recommended dilution:	500-2000

