

小鼠抗 NQO2 单克隆抗体

中文名称: 小鼠抗 NQO2 单克隆抗体

英文名称: Anti-NQO2 mouse monoclonal antibody

别名: DHQV; DIA6; NMOR2; QR2

储存: 冷冻 (-20℃) 避光

抗原: NQO2

宿主: Mouse

反应种属: Human , Dog , Monkey

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	NQO2 (EC 1.10.99.2) is a flavoprotein that catalyzes the 2-electron reduction of various quinones, redox dyes, and the vitamin K menadiolone. NQO2 predominantly uses dihydronicotinamide riboside (NRH) as the electron donor (summary by Wu et al., 1997 [PubMed 9367528]). [supplied by OMIM].
Applications:	WB, IHC, IF
Name of antibody:	NQO2
Immunogen:	Fusion protein of human NQO2
Full name:	NAD(P)H dehydrogenase, quinone 2 (NQO2), transcript variant 3
Synonyms:	DHQV; DIA6; NMOR2; QR2
SwissProt:	P16083

IHC positive control:	human pancreas tissue and human liver tissue; carcinoma of human liver tissue
IHC Recommended dilution:	30-150
WB Predicted band size:	26 kDa
WB Positive control:	HepG2, HeLa, SVT2, A549, COS7, Jurkat, MDCK, MCF-7 cell lysates
WB Recommended dilution:	500-2000