

兔抗 TENT5C 多克隆抗体

中文名称：兔抗 TENT5C 多克隆抗体

英文名称：Anti-TENT5C rabbit polyclonal antibody

别名：terminal nucleotidyltransferase 5C; FAM46C

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TENT5C

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Probable nucleotidyltransferase that may act as a non-canonical poly(A) RNA polymerase. (Microbial infection) Seems to enhance replication of some viruses, including yellow fever virus, in response to type I interferon.
Applications:	ELISA, IHC
Name of antibody:	TENT5C
Immunogen:	Fusion protein of human TENT5C
Full name:	terminal nucleotidyltransferase 5C
Synonyms:	FAM46C
SwissProt:	Q5VWP2
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	150-300

