

小鼠抗 PAAF1 单克隆抗体

中文名称: 小鼠抗 PAAF1 单克隆抗体

英文名称: Anti-PAAF1 mouse monoclonal antibody

别名: PAAF; Rpn14; WDR71

抗原: PAAF1

储存: 冷冻 (-20°C) 避光

宿主: Mouse

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes a WD repeat-containing protein involved in regulation of association of proteasome components. During HIV infection, the encoded protein is thought to promote provirus transcription through recruitment of the 19S regulatory complex. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	WB, IF
Name of antibody:	PAAF1
Immunogen:	Fusion protein of human PAAF1
Full name:	proteasomal ATPase-associated factor 1 (PAAF1), transcript variant 4

Synonyms:	PAAF; Rpn14; WDR71
SwissProt:	Q9BRP4
WB Predicted band size :	42 kDa
WB Positive control:	HepG2, Hela, SVT2 cell lysates
WB Recommended dilution:	200-1000