

小鼠抗 HDAC6 单克隆抗体

中文名称： 小鼠抗 HDAC6 单克隆抗体

英文名称： Anti-HDAC6 mouse monoclonal antibody

别名： histone deacetylase 6; HD6; JM21; CPBHM; PPP1R90

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： HDAC6

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene belongs to class II of the histone deacetylase/acuc/apha family. It contains an internal duplication of two catalytic domains which appear to function independently of each other. This protein possesses histone deacetylase activity and represses transcription.
Applications:	WB
Name of antibody:	HDAC6

Immunogen:	Fusion protein of human HDAC6
Full name:	histone deacetylase 6
Synonyms:	HD6; JM21; CPBHM; PPP1R90
SwissProt:	Q9UBN7
WB Predicted band size:	131 kDa
WB Positive control:	T47D cell, Jurkat cell, CHO-K1 cell and CHO-K1 transfected by HDAC6-fragment EGFP fusion protein cell lysates
WB Recommended dilution:	1000-5000

