

小鼠抗 Pan-CK 单克隆抗体

中文名称： 小鼠抗 Pan-CK 单克隆抗体

英文名称： Anti-Pan-CK mouse monoclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃） 避光

宿 主： Mouse

抗 原： Pan-CK

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	Cytokeratins, a group comprising at least 29 different proteins, are characteristic of epithelial and trichocytic cells. Cytokeratins 1, 4, 5, 6, and 8 are members of the type II neutral to basic subfamily. Monoclonal anti cytokeratins are specific markers of epithelial cell differentiation and have been widely used as tools in tumor identification and classification.
Applications:	IHC
Name of antibody:	Pan-CK
Immunogen:	Human epidermal cytokeratins
Full name:	pan-cytokeratin
IHC positive control:	adenocarcinoma of human breast tissue and human colon tissue and adenocarcinoma of human colon tissue and human Kidney tissue

IHC Recommend dilution:

30-150