

小鼠抗 HDHD1 单克隆抗体

中文名称： 小鼠抗 HDHD1 单克隆抗体

英文名称： Anti-HDHD1 mouse monoclonal antibody

别名： GS1; HDHD1; HDHD1A; FAM16AX; DXF68S1E

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： HDHD1

反应种属： Human, Dog, Rat, Monkey

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a member of the haloacid dehalogenase-like (HAD) hydrolase superfamily. The encoded protein has no known biological function. This gene has a pseudogene on chromosome 1. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq].
Applications:	WB, IHC, FC
Name of antibody:	HDHD1
Immunogen:	Fusion protein of human HDHD1
Full name:	haloacid dehalogenase-like hydrolase domain containing 1 (HDHD1), transcript variant 2

Synonyms:	GS1; HDHD1; HDHD1A; FAM16AX; DXF68S1E
SwissProt:	Q08623
IHC positive control:	adenocarcinoma of human endometrium tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	25 kDa
WB Positive control:	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
WB Recommended dilution:	500-2000