

小鼠抗 PARVA 单克隆抗体

中文名称: 小鼠抗 PARVA 单克隆抗体

英文名称: Anti-PARVA mouse monoclonal antibody

别名: CH-ILKBP; MXRA2

抗原: PARVA

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human , Dog , Rat , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes a member of the parvin family of actin-binding proteins. Parvins are associated with focal contacts and contain calponin homology domains that bind to actin filaments. The encoded protein is part of the integrin-linked kinase signaling complex and plays a role in cell adhesion, motility and survival.
Applications:	WB, IHC, IF
Name of antibody:	PARVA
Immunogen:	Fusion protein of human PARVA
Full name:	parvin, alpha (PARVA)
Synonyms:	CH-ILKBP; MXRA2

SwissProt:	Q9NVD7
IHC positive control:	adenocarcinoma of human breast tissue and human kidney tissue; human pancreas tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	42 kDa
WB Positive control:	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
WB Recommended dilution:	200-1000