

小鼠抗 HOMER1 单克隆抗体

中文名称：小鼠抗 HOMER1 单克隆抗体

英文名称：Anti-HOMER1 mouse monoclonal antibody

别名：homer scaffolding protein 1; HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Mouse

抗原：HOMER1

反应种属：Human, Rat, Mouse

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	This gene encodes a member of the homer family of dendritic proteins. Members of this family regulate group 1 metabotropic glutamate receptor function.
Applications:	WB, IHC
Name of antibody:	HOMER1
Immunogen:	Fusion protein of human HOMER1
Full name:	homer scaffolding protein 1
Synonyms:	HOMER; SYN47; Ves-1; HOMER1A; HOMER1B; HOMER1C
SwissProt:	Q86YM7
IHC positive control:	Human spleen tissue

IHC Recommend dilution:	200-500
WB Predicted band size:	40 KD
WB Positive control:	HepG2 cell, Mouse brain tissue, Rat brain tissue lysates
WB Recommended dilution:	500-2000