

小鼠抗 HOXA3 单克隆抗体

中文名称：小鼠抗 HOXA3 单克隆抗体

英文名称：Anti-HOXA3 mouse monoclonal antibody

别名：homeobox A3; HOX1; HOX1E

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：HOXA3

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation. Three transcript variants encoding two different isoforms have been found for this gene.
Applications:	WB

Name of antibody:	HOXA3
Immunogen:	Fusion protein of human HOXA3
Full name:	homeobox A3
Synonyms:	HOX1; HOX1E
SwissProt:	O43365
WB Predicted band size:	46 KD
WB Positive control:	HepG2 cell lysate
WB Recommended dilution:	500-2000