

小鼠抗 PBLD 单克隆抗体

- 中文名称： 小鼠抗 PBLD 单克隆抗体
- 英文名称： Anti-PBLD mouse monoclonal antibody
- 相关类别： 一抗
- 储 存： 冷冻（-20℃）
- 宿 主： Mouse
- 抗 原： PBLD
- 反应种属： Human
- 标 记 物： Unconjugate
- 克隆类型： mouse monoclonal

技术规格

Background:	MAWDBP (MAWD binding protein), also known as PBLD (phenazine biosynthesis-like protein domain containing) or MAWBP, is a 288 amino acid protein that belongs to the phenazine biosynthesis-like protein (PhzF) family. It has been suggested that MAWDBP is the only representative of the PhzF family in the human genome. Expressed in most tissues, MAWDBP is a WD-40 repeat-containing 2-propeller protein believed to participate in the MAPK signaling pathway. Involved in multiple basic cellular functions, expression of MAWDBP is elevated in several disease processes, including insulin resistance, folate deficiency and hypotension. It is thought that MAWDBP may also be involved in carcinogenesis.
Applications:	ELISA, WB

Name of antibody:	PBLD
Immunogen:	Fusion protein of human PBLD
Full name:	phenazine biosynthesis like protein domain containing
Synonyms:	MAWBP; MAWDBP; HEL-S-306; MAWBP(PBLD)
SwissProt:	P30039
WB Predicted band size :	32 kDa
WB Positive control:	BGC823 whole cell lysate
WB Recommended dilution:	1000-5000

