

## 小鼠抗 HS1BP3 单克隆抗体

中文名称：小鼠抗 HS1BP3 单克隆抗体

英文名称：Anti-HS1BP3 mouse monoclonal antibody

别名：HCLS1 binding protein 3; ETM2; HS1-BP3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：HS1BP3

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene shares similarity with mouse Hs1bp3, an Hcls1/Hs1-interacting protein that may be involved in lymphocyte activation.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	HS1BP3
<b>Immunogen:</b>	Fusion protein of human HS1BP3
<b>Full name:</b>	HCLS1 binding protein 3
<b>Synonyms:</b>	ETM2; HS1-BP3
<b>SwissProt:</b>	Q53T59
<b>IHC positive control:</b>	Human kidney tissue
<b>IHC Recommend dilution:</b>	200-500
<b>WB Predicted band size:</b>	43 KD



<b>WB Positive control:</b>	Hela, 293T, Jurkat cell lysates
<b>WB Recommended dilution:</b>	500-2000