

## 小鼠抗 PCDH7 单克隆抗体

中文名称: 小鼠抗 PCDH7 单克隆抗体

英文名称: Anti-PCDH7 mouse monoclonal antibody

别名: BH-Pcdh; BHPCDH; PPP1R120

抗原: PCDH7

储存: 冷冻 (-20°C) 避光

宿主: Mouse

反应种属: Human , Dog , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	This gene belongs to the protocadherin gene family, a subfamily of the cadherin superfamily. The gene encodes a protein with an extracellular domain containing 7 cadherin repeats. The gene product is an integral membrane protein that is thought to function in cell-cell recognition and adhesion. Alternative splicing yields isoforms with unique cytoplasmic tails. [provided by RefSeq, Jul 2008].
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	PCDH7
<b>Immunogen:</b>	Fusion protein of human PCDH7
<b>Full name:</b>	protocadherin 7 (PCDH7), transcript variant c

<b>Synonyms:</b>	BH-Pcdh; BHPCDH; PPP1R120
<b>SwissProt:</b>	O60245
<b>IHC positive control:</b>	human endometrium tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	116 kDa
<b>WB Positive control:</b>	HepG2, Hela, SVT2, COS7, MDCK, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	200-1000