

## 小鼠抗 HSF1 单克隆抗体

中文名称： 小鼠抗 HSF1 单克隆抗体

英文名称： Anti-HSF1 mouse monoclonal antibody

相关类别： 一抗

储 存： 冷冻 (-20℃)

宿 主： Mouse

抗 原： HSF1

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	The product of this gene is a transcription factor that is rapidly induced after temperature stress and binds heat shock promoter elements (HSE). This protein plays a role in the regulation of lifespan. Expression of this gene is repressed by phosphorylation, which promotes binding by heat shock protein 90.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	HSF1
<b>Immunogen:</b>	Fusion protein of human HSF1
<b>Full name:</b>	heat shock transcription factor 1
<b>Synonyms:</b>	HSTF1
<b>SwissProt:</b>	Q00613
<b>WB Predicted band size:</b>	57 kDa

<b>WB Positive control:</b>	HeLa whole cell lysate
<b>WB Recommended dilution:</b>	500-2000

