

## 兔抗 TGM3 多克隆抗体

中文名称：兔抗 TGM3 多克隆抗体

英文名称： Anti-TGM3 rabbit polyclonal antibody

别名： transglutaminase 3; TGE; UHS2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TGM3

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Transglutaminases are enzymes that catalyze the crosslinking of proteins by epsilon-gamma glutamyl lysine isopeptide bonds. While the primary structure of transglutaminases is not conserved, they all have the same amino acid sequence at their active sites and their activity is calcium-dependent. The protein encoded by this gene consists of two polypeptide chains activated from a single precursor protein by proteolysis. The encoded protein is involved in the later stages of cell envelope formation in the epidermis and hair follicle.
<b>Applications:</b>	ELISA, WB, IHC

<b>Name of antibody:</b>	TGM3
<b>Immunogen:</b>	Fusion protein of human TGM3
<b>Full name:</b>	transglutaminase 3
<b>Synonyms:</b>	TGE; UHS2
<b>SwissProt:</b>	Q08188
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	50- 200
<b>WB Predicted band size:</b>	77 kDa
<b>WB Positive control:</b>	Human tongue tissue lysate
<b>WB Recommended dilution:</b>	1000-5000



