

## 兔抗 TGM5 多克隆抗体

中文名称：兔抗 TGM5 多克隆抗体

英文名称： Anti-TGM5 rabbit polyclonal antibody

别名： TGX; PSS2; TGM6; TGMX; TGASE5; TGASEX

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TGM5

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the transglutaminase family. The encoded protein catalyzes formation of protein cross-links between glutamine and lysine residues, often resulting in stabilization of protein assemblies. This reaction is calcium dependent. Mutations in this gene have been associated with acral peeling skin syndrome.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	TGM5
<b>Immunogen:</b>	Fusion protein of human TGM5
<b>Full name:</b>	transglutaminase 5
<b>Synonyms :</b>	TGX; PSS2; TGM6; TGMX; TGASE5; TGASEX

<b>SwissProt:</b>	O43548
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and human ovarian cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	81 kDa
<b>WB Positive control:</b>	Hela and A431 cells, Mouse skin tissue, MCF7 and K562 cells
<b>WB Recommended dilution:</b>	500-2000



