

兔抗 TGM6 多克隆抗体

中文名称：兔抗 TGM6 多克隆抗体

英文名称：Anti-TGM6 rabbit polyclonal antibody

别名：TG6; TGY; SCA35; TGM3L; dJ734P14.3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TGM6

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the transglutaminase superfamily. It catalyzes the cross-linking of proteins and the conjugation of polyamines to proteins. Mutations in this gene are associated with spinocerebellar ataxia type 35 (SCA35). Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	TGM6
Immunogen:	Synthetic peptide of human TGM6
Full name:	transglutaminase 6

Synonyms :	TG6; TGY; SCA35; TGM3L; dJ734P14.3
SwissProt:	O95932
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	25-100

