

小鼠抗 PECAM1 单克隆抗体

- 中文名称: 小鼠抗 PECAM1 单克隆抗体
- 英文名称: Anti-PECAM1 mouse monoclonal antibody
- 别名: CD31; CD31/EndoCAM; endoCAM; GPIIA ; PECA1; PECAM-1
- 储存: 冷冻 (-20℃) 避光
- 抗原: PECAM1
- 宿主: Mouse
- 反应种属: Human
- 标记物: Unconjugate
- 相关类别: 一抗
- 克隆类型: mouse monoclonal

技术规格

Background:	The protein encoded by this gene is found on the surface of platelets, monocytes, neutrophils, and some types of T-cells, and makes up a large portion of endothelial cell intercellular junctions. The encoded protein is a member of the immunoglobulin superfamily and is likely involved in leukocyte migration, angiogenesis, and integrin activation. [provided by RefSeq].
Applications:	WB, IHC, IF, FC
Name of antibody:	PECAM1
Immunogen:	Fusion protein of human PECAM1
Full name:	platelet/endothelial cell adhesion molecule 1 (PECAM1)
Synonyms:	CD31; CD31/EndoCAM; endoCAM; GPIIA ; PECA1; PECAM-1

SwissProt:	P16284
IHC positive control:	human pancreas tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	83 kDa
WB Positive control:	Jurkat and PC3 cell lysates
WB Recommended dilution:	500-2000



