

小鼠抗 PGAM1 单克隆抗体

中文名称： 小鼠抗 PGAM1 单克隆抗体

英文名称： Anti-PGAM1 mouse monoclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Mouse

抗 原： PGAM1

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a mutase that catalyzes the reversible reaction of 3-phosphoglycerate (3-PGA) to 2-phosphoglycerate (2-PGA) in the glycolytic pathway. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	PGAM1
Immunogen:	Fusion protein of human PGAM1
Full name:	phosphoglycerate mutase 1
Synonyms:	PGAMA; PGAM-B; HEL-S-35
SwissProt:	P18669
WB Predicted band size:	29 kDa
WB Positive control:	Jurkat whole cell lysate

WB Recommended dilution:	500-2000
--------------------------	----------

