

## 兔抗 THAP2 多克隆抗体

中文名称：兔抗 THAP2 多克隆抗体

英文名称： Anti-THAP2 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： THAP2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	THAP2, also named as THAP domain-containing protein 2, is a 228 amino acid protein, which contains one THAP-type zinc finger. Members of the THAP (thanatos-associated protein) family of proteins contain a well conserved DNA-binding domain known as the THAP-type zinc finger motif. Proteins containing the THAP-type zinc finger motif are commonly involved in transcriptional regulation, cell cycle control, apoptosis and chromatin modification
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	THAP2
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	THAP domain containing, apoptosis associated prot

	ein 2
<b>SwissProt:</b>	Q9H0W7
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

