

## 小鼠抗 HSPA8 单克隆抗体

中文名称： 小鼠抗 HSPA8 单克隆抗体

英文名称： Anti-HSPA8 mouse monoclonal antibody

别名： heat shock protein family A (Hsp70) member 8; LAP1; HSC54; HSC70; HSC71; HSP71; HSP73; LAP-1; NIP71; HEL-33; HSPA10; HEL-S-72p

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： HSPA8

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

**Background:**

This gene encodes a member of the heat shock protein 70 family, which contains both heat-inducible and constitutively expressed members. This protein belongs to the I<sub>α</sub> group, which are also referred to as heat-shock cognate proteins. It functions as a chaperone, and binds to nascent polypeptides to facilitate correct folding. It also functions as an ATPase in the disassembly of clathrin-coated vesicles during transport of membrane components through the cell. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

	ne.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	HSPA8
<b>Immunogen:</b>	Fusion protein of Human HSPA8
<b>Full name:</b>	heat shock protein family A (Hsp70) member 8
<b>Synonyms:</b>	LAP1; HSC54; HSC70; HSC71; HSP71; HSP73; LAP-1; NIP7 1; HEL-33; HSPA10; HEL-S-72p
<b>SwissProt:</b>	P11142
<b>WB Predicted band size:</b>	71 kDa
<b>WB Positive control:</b>	A549 and Hela cell, Rat brain tissue, Mouse brain tissue, 293T, Raji and NIH/3T3 cell lysates
<b>WB Recommended dilution:</b>	2000-10000

