

兔抗 THAP10 多克隆抗体

中文名称：兔抗 THAP10 多克隆抗体

英文名称：Anti-THAP10 rabbit polyclonal antibody

抗 原：THAP10

宿 主：Rabbit

反应种属：Human

相关类别：一抗

储 存：冷冻（-20℃）

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of a family of proteins sharing an N-terminal Thanatos-associated domain. The Thanatos-associated domain contains a zinc finger signature similar to DNA-binding domains. This gene is part of a bidirectional gene pair on the long arm of chromosome 15 that is regulated by estrogen and may play a role in breast cancer.
Applications:	ELISA, IHC
Name of antibody:	THAP10
Immunogen:	Full length fusion protein

Full name:	THAP domain containing 10
SwissProt:	Q9P2Z0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	30-150

