

小鼠抗 ICAM1 单克隆抗体

- 中文名称: 小鼠抗 ICAM1 单克隆抗体
- 英文名称: Anti-ICAM1 mouse monoclonal antibody
- 别名: intercellular adhesion molecule 1; BB2; CD54; P3.58
- 相关类别: 一抗
- 储存: 冷冻 (-20℃)
- 宿主: Mouse
- 抗原: ICAM1
- 反应种属: Human
- 标记物: Unconjugate
- 克隆类型: Mouse Monoclonal

技术规格

Background:	This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.
Applications:	WB, IF
Name of antibody:	ICAM1
Immunogen:	Fusion protein of human ICAM1
Full name:	intercellular adhesion molecule 1
Synonyms:	BB2; CD54; P3.58
SwissProt:	P05362

WB Predicted band size:	58 kDa
WB Positive control:	SW480 and Raji cell lysates
WB Recommended dilution:	1000-5000
IF positive control:	SW480 cells
IF Recommend dilution:	100-200



