

## 兔抗 THBS1 多克隆抗体

中文名称：兔抗 THBS1 多克隆抗体

英文名称：Anti-THBS1 rabbit polyclonal antibody

别名：TSP; THBS; TSP1; TSP-1; THBS-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：THBS1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a subunit of a disulfide-linked homotrimeric protein. This protein is an adhesive glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. This protein can bind to fibrinogen, fibronectin, laminin, type V collagen and integrins alpha-V/beta-1. This protein has been shown to play roles in platelet aggregation, angiogenesis, and tumorigenesis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	THBS1
<b>Immunogen:</b>	Synthetic peptide of human THBS1

<b>Full name:</b>	thrombospondin 1
<b>Synonyms :</b>	TSP; THBS; TSP1; TSP-1; THBS-1
<b>SwissProt:</b>	P07996
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human cervical cancer
<b>IHC Recommend dilution:</b>	50-200

