

## 兔抗 THBS2 多克隆抗体

中文名称：兔抗 THBS2 多克隆抗体

英文名称：Anti-THBS2 rabbit polyclonal antibody

别名：TSP2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：THBS2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the thrombospondin family. It is a disulfide-linked homotrimeric glycoprotein that mediates cell-to-cell and cell-to-matrix interactions. This protein has been shown to function as a potent inhibitor of tumor growth and angiogenesis. Studies of the mouse counterpart suggest that this protein may modulate the cell surface properties of mesenchymal cells and be involved in cell adhesion and migration.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	THBS2

<b>Immunogen:</b>	Synthetic peptide of humanTHBS2
<b>Full name:</b>	thrombospondin 2
<b>Synonyms :</b>	TSP2
<b>SwissProt:</b>	P35442
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

