

小鼠抗 PIK3R2 单克隆抗体

中文名称： 小鼠抗 PIK3R2 单克隆抗体

英文名称： Anti-PIK3R2 mouse monoclonal antibody

别名： phosphoinositide-3-kinase regulatory subunit 2; p85; MPPH; P85B; MPPH1; p85-BETA

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： PIK3R2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	Phosphatidylinositol 3-kinase (PI3K) is a lipid kinase that phosphorylates phosphatidylinositol and similar compounds, creating second messengers important in growth signaling pathways. PI3K functions as a heterodimer of a regulatory and a catalytic subunit. The protein encoded by this gene is a regulatory component of PI3K. Two transcript variants, one protein coding and the other non-protein coding, have been found for this gene.
Applications:	WB
Name of antibody:	PIK3R2

Immunogen:	Fusion protein of human PIK3R2
Full name:	phosphoinositide-3-kinase regulatory subunit 2
Synonyms:	p85; MPPH; P85B; MPPH1; p85-BETA
SwissProt:	O00459
WB Predicted band size:	82 kDa
WB Positive control:	C6, K562, MCF7, U87MG, HCT116, 3T3 and HeLa cell lysates
WB Recommended dilution :	1000-5000

