

## 兔抗 THOP1 多克隆抗体

中文名称：兔抗 THOP1 多克隆抗体

英文名称：Anti-THOP1 rabbit polyclonal antibody

别名：thimet oligopeptidase 1; TOP; MP78; EP24.15; MEPD\_HUMAN

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：THOP1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a kinase that uses zinc as a cofactor. The encoded oligopeptidase cleaves cytosolic peptides, making them unavailable for display on antigen-presenting cells. This protein also cleaves neuropeptides under 20 aa in length and can degrade beta-amyloid precursor protein to amyloidogenic peptides.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	THOP1
<b>Immunogen:</b>	Fusion protein of human THOP1
<b>Full name:</b>	thimet oligopeptidase 1

<b>Synonyms:</b>	TOP; MP78; EP24.15; MEPD_HUMAN
<b>SwissProt:</b>	P52888
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	79 kDa
<b>WB Positive control:</b>	Human cerebrum tissue and Human fetal brain tissue lysates
<b>WB Recommended dilution:</b>	500-2000



