

小鼠抗 PLAC8 单克隆抗体

中文名称： 小鼠抗 PLAC8 单克隆抗体

英文名称： Anti-PLAC8 mouse monoclonal antibody

别名： placenta specific 8; C15; DGIC; onzin; PNAS-144

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： PLAC8

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	PLAC8, also named as Protein C15, is a placenta-specific protein. The expression was downregulated upon activation, particularly in dendritic cells generated in vitro from CD34 (142 230) progenitor cells. The MW of PLAC8 is about 12-20kd
Applications:	WB, IHC
Name of antibody:	PLAC8
Immunogen:	Fusion protein of human PLAC8
Full name:	placenta specific 8
Synonyms:	C15; DGIC; onzin; PNAS-144
SwissProt:	Q9NZF1
IHC positive control:	Human tonsil tissue

IHC Recommend dilution:	200-500
WB Predicted band size:	13 KD
WB Positive control:	Human liver tissue lysate
WB Recommended dilution :	200-1000