

小鼠抗 IQGAP1 单克隆抗体

中文名称： 小鼠抗 IQGAP1 单克隆抗体

英文名称： Anti-IQGAP1 mouse monoclonal antibody

别名： IQ motif containing GTPase activating protein 1; SAR1; p195; HUMORFA01

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： IQGAP1

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a member of the IQGAP family. The protein contains four IQ domains, one calponin homology domain, one Ras-GAP domain and one WW domain. It interacts with components of the cytoskeleton, with cell adhesion molecules, and with several signaling molecules to regulate cell morphology and motility. Expression of the protein is upregulated by gene amplification in two gastric cancer cell lines.
Applications:	WB, IHC
Name of antibody:	IQGAP1
Immunogen:	Fusion protein of human IQGAP1

Full name:	IQ motif containing GTPase activating protein 1
Synonyms:	SAR1; p195; HUMORFA01
SwissProt:	P46940
IHC positive control:	Human colorectal cancer and Human pancreas tissue
IHC Recommend dilution:	500-2000
WB Predicted band size:	189 KD
WB Positive control:	Hela, 293T, Jurkat, MCF7, K562, A549, HepG2 cell lysates
WB Recommended dilution:	500-2000