

小鼠抗 PMS1 单克隆抗体

中文名称： 小鼠抗 PMS1 单克隆抗体

英文名称： Anti-PMS1 mouse monoclonal antibody

别名： PMS1 homolog 1, mismatch repair system component; MLH2; PMSL1; hPMS1; HNPCC3

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： PMS1

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a protein belonging to the DNA mismatch repair mutL/hexB family. This protein is thought to be involved in the repair of DNA mismatches, and it can form heterodimers with MLH1, a known DNA mismatch repair protein. Mutations in this gene cause hereditary nonpolyposis colorectal cancer type 3 (HNPCC3) either alone or in combination with mutations in other genes involved in the HNPCC phenotype, which is also known as Lynch syndrome.
Applications:	WB
Name of antibody:	PMS1

Immunogen:	Fusion protein of human PMS1
Full name:	PMS1 homolog 1, mismatch repair system component
Synonyms:	MLH2; PMSL1; hPMS1; HNPCC3
SwissProt:	P54277
WB Predicted band size:	106 KD
WB Positive control:	Hela, Jurkat, HepG2 cell lysates
WB Recommended dilution:	200-1000