

小鼠抗 IST1 单克隆抗体

中文名称： 小鼠抗 IST1 单克隆抗体

英文名称： Anti-IST1 mouse monoclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Mouse

抗 原： IST1

反应种属： Human

标 记 物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a protein with MIT-interacting motifs that interacts with components of endosomal sorting complexes required for transport (ESCRT). ESCRT functions in vesicle budding, such as that which occurs during membrane abscission in cytokinesis. There is a pseudogene for this gene on chromosome 19. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, WB
Name of antibody:	IST1
Immunogen:	Fusion protein of human IST1
Full name:	IST1, ESCRT-III associated factor
Synonyms:	OLC1; hIST1

SwissProt:	P53990
WB Predicted band size:	40 kDa
WB Positive control:	H520 whole cell lysate
WB Recommended dilution:	500-2000

