

小鼠抗 POGK 单克隆抗体

中文名称: 小鼠抗 POGK 单克隆抗体

英文名称: Anti-POGK mouse monoclonal antibody

别名: BASS2; KRBOX2; LST003

储存: 冷冻 (-20℃) 避光

抗原: POGK

宿主: Mouse

反应种属: Human , Monkey

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	The exact function of the protein encoded by this gene is not known. However, this gene product contains a KRAB domain (which is involved in protein-protein interactions) at the N-terminus, and a transposase domain at the C-terminus, suggesting that it may belong to the family of DNA-mediated transposons in human. [provided by RefSeq, Jul 2008].
Applications:	WB, IHC, IF
Name of antibody:	POGK
Immunogen:	Fusion protein of human POGK
Full name:	pogo transposable element with KRAB domain (POGK)
Synonyms:	BASS2; KRBOX2; LST003

SwissProt:	Q9P215
IHC positive control:	human kidney tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	69 kDa
WB Positive control:	HepG2, Hela, A549, COS7, Jurkat, MCF-7 cell lysates
WB Recommended dilution:	500-2000