

兔抗 THUMPD3 多克隆抗体

中文名称：兔抗 THUMPD3 多克隆抗体

英文名称：Anti-THUMPD3 rabbit polyclonal antibody

别名：THUMP domain containing 3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：THUMPD3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	THUMPD3
Immunogen:	Fusion protein of human THUMPD3
Full name:	THUMP domain containing 3
SwissProt:	Q9BV44
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human prostate cancer and Human colorectal cancer
IHC Recommend dilution:	100-300

