

小鼠抗 ITGB3(N-term)单克隆抗体

中文名称： 小鼠抗 ITGB3(N-term)单克隆抗体

英文名称： Anti-ITGB3(N-term) mouse monoclonal antibody

别名： integrin subunit beta 3; GT; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： ITGB3(N-term)

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	The ITGB3 protein product is the integrin beta chain beta 3. Integrins are integral cell-surface proteins composed of an alpha chain and a beta chain. A given chain may combine with multiple partners resulting in different integrins. Integrin beta 3 is found along with the alpha IIb chain in platelets. Integrins are known to participate in cell adhesion as well as cell-surface mediated signalling.
Applications:	WB
Name of antibody:	ITGB3(N-term)
Immunogen:	Fusion protein of human ITGB3
Full name:	integrin subunit beta 3

Synonyms:	GT; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16
SwissProt:	P05106
WB Predicted band size:	87 kDa
WB Positive control:	15ug and 25ug Leukocyte whole cell lysates
WB Recommended dilution:	1000-5000

