

## 小鼠抗 POLR2E 单克隆抗体

中文名称： 小鼠抗 POLR2E 单克隆抗体

英文名称： Anti-POLR2E mouse monoclonal antibody

别名： hRPB25; hsRPB5; RPABC1; RPB5; XAP4

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： POLR2E

反应种属： Human , Dog , Rat , Monkey , Mouse

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes the fifth largest subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. This subunit is shared by the other two DNA-directed RNA polymerases and is present in two-fold molar excess over the other polymerase subunits. An interaction between this subunit and a hepatitis virus transactivating protein has been demonstrated, suggesting that interaction between transcriptional activators and the polymerase can occur through this subunit. A pseudogene is located on chromosome 11. [provided by RefSeq].
<b>Applications:</b>	WB, IHC, IF, FC

<b>Name of antibody:</b>	POLR2E
<b>Immunogen:</b>	Fusion protein of human POLR2E
<b>Full name:</b>	polymerase (RNA) II (DNA directed) polypeptide E, 25kDa (POLR2E)
<b>Synonyms:</b>	hRPB25; hsRPB5; RPABC1; RPB5; XAP4
<b>SwissProt:</b>	P19388
<b>IHC positive control:</b>	human breast tissue and carcinoma of human lung tissue; human pancreas tissue and carcinoma of human prostate tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	25 kDa
<b>WB Positive control:</b>	HepG2, HeLa, SVT2, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	500-2000