

小鼠抗 ITM2A 单克隆抗体

中文名称：小鼠抗 ITM2A 单克隆抗体

英文名称：Anti-ITM2A mouse monoclonal antibody

别名：integral membrane protein 2A; E25A; BRICD2A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：ITM2A

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	This gene encodes a type II membrane protein that belongs to the ITM2 family. Studies in mouse suggest that it may be involved in osteo- and chondrogenic differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	ITM2A
Immunogen:	Fusion protein of human ITM2A
Full name:	integral membrane protein 2A
Synonyms:	E25A; BRICD2A
SwissProt:	O43736



WB Predicted band size:	30 KD
WB Positive control:	K562 cell lysate
WB Recommended dilution :	500-2000