

小鼠抗 PPIF 单克隆抗体

- 中文名称： 小鼠抗 PPIF 单克隆抗体
- 英文名称： Anti-PPIF mouse monoclonal antibody
- 别名： peptidylprolyl isomerase F; CYP3; CypD; CyP-M; Cyp-D
- 相关类别： 一抗
- 储存： 冷冻（-20℃）
- 宿主： Mouse
- 抗原： PPIF
- 反应种属： Human
- 标记物： Unconjugate
- 克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a member of the peptidyl-prolyl cis-trans isomerase (PPIase) family. PPIases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and accelerate the folding of proteins. This protein is part of the mitochondrial permeability transition pore in the inner mitochondrial membrane. Activation of this pore is thought to be involved in the induction of apoptotic and necrotic cell death.
Applications:	WB, IHC
Name of antibody:	PPIF
Immunogen:	Fusion protein of human PPIF
Full name:	peptidylprolyl isomerase F

Synonyms:	CYP3; CypD; CyP-M; Cyp-D
SwissProt:	P30405
IHC positive control:	Human liver cancer
IHC Recommend dilution:	500-1000
WB Predicted band size:	22 KD
WB Positive control:	A431, Hela, HepG2 and K562 cell lysates
WB Recommended dilution:	500-2000