

兔抗 TIGD5 多克隆抗体

中文名称：兔抗 TIGD5 多克隆抗体

英文名称： Anti-TIGD5 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： TIGD5

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. The exact function of this gene is not known.
Applications:	ELISA, IHC
Name of antibody:	TIGD5
Immunogen:	Fusion protein of human TIGD5
Full name:	tigger transposable element derived 5
SwissProt:	Q53EQ6
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	25-100



