

小鼠抗 KRT5 单克隆抗体

- 中文名称： 小鼠抗 KRT5 单克隆抗体
- 英文名称： Anti-KRT5 mouse monoclonal antibody
- 别名： keratin 5; K5; CK5; DDD; DDD1; EBS2; KRT5A
- 相关类别： 一抗
- 储存： 冷冻（-20℃）
- 宿主： Mouse
- 抗原： KRT5
- 反应种属： Human, Mouse, Rat
- 标记物： Unconjugate
- 克隆类型： Mouse Monoclonal

技术规格

Background:

The protein encoded by this gene is a member of the keratin gene family. The type II cytokeratins consist of basic or neutral proteins which are arranged in pairs of heterotypic keratin chains coexpressed during differentiation of simple and stratified epithelial tissues. This type II cytokeratin is specifically expressed in the basal layer of the epidermis with family member KRT14. Mutations in these genes have been associated with a complex of diseases termed epidermolysis bullosa simplex. The type II cytokeratins are clustered in a region of chromosome 12q12-q13.

Applications:	WB
Name of antibody:	KRT5
Immunogen:	Fusion protein of human KRT5
Full name:	keratin 5
Synonyms:	K5; CK5; DDD; DDD1; EBS2; KRT5A
SwissProt:	P13647
WB Predicted band size:	62 kDa
WB Positive control:	C6, 3T3, A549, MCF7 and Hela cell lysates
WB Recommended dilution:	5000-10000

