

## 小鼠抗 PPP2CA/PPP2CB 单克隆抗体

中文名称： 小鼠抗 PPP2CA/PPP2CB 单克隆抗体

英文名称： Anti-PPP2CA/PPP2CB mouse monoclonal antibody

别名： protein phosphatase 2 catalytic subunit alpha/beta; RP-C; PP2Ac; PP2CA; PP2Calpha; PP2CB; PP2Abeta

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： PPP2CA/PPP2CB

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： Mouse Monoclonal

### 技术规格

<b>Background:</b>	This gene encodes the phosphatase 2A catalytic subunit. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. This gene encodes an alpha isoform and a beta isoform of the catalytic subunit.
<b>Applications:</b>	WB, IP
<b>Name of antibody:</b>	PPP2CA/PPP2CB
<b>Immunogen:</b>	Fusion protein of human PPP2CA/PPP2CB

<b>Full name:</b>	protein phosphatase 2 catalytic subunit alpha/beta
<b>Synonyms:</b>	RP-C; PP2Ac; PP2CA; PP2Calpha; PP2CB; PP2Abeta
<b>SwissProt:</b>	P67775/P62714
<b>WB Predicted band size:</b>	36 kDa
<b>WB Positive control:</b>	3T3, C6, MCF7 and Hela cell lysates
<b>WB Recommended dilution:</b>	1000-5000



