

兔抗 TIMP4 多克隆抗体

中文名称：兔抗 TIMP4 多克隆抗体

英文名称： Anti-TIMP4 rabbit polyclonal antibody

抗 原： TIMP4

储 存： 冷冻（-20℃）

宿 主： Rabbit

反应种属： Human, Mouse, Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene belongs to the TIMP gene family. The proteins encoded by this gene family are inhibitors of the matrix metalloproteinases, a group of peptidases involved in degradation of the extracellular matrix. The secreted, netrin domain-containing protein encoded by this gene is involved in regulation of platelet aggregation and recruitment and may play role in hormonal regulation and endometrial tissue remodeling.
Applications:	ELISA, WB, IHC
Name of antibody:	TIMP4
Immunogen:	Fusion protein of human TIMP4
Full name:	TIMP metalloproteinase inhibitor 4

SwissProt:	Q99727
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	15-50
WB Predicted band size:	26 kDa
WB Positive control:	Human brain tissue, PC3 and hela cells
WB Recommended dilution:	200-1000



