

小鼠抗 KRT8 单克隆抗体

中文名称： 小鼠抗 KRT8 单克隆抗体

英文名称： Anti-KRT8 mouse monoclonal antibody

别名： keratin 8; K8; KO; CK8; CK-8; CYK8; K2C8; CARD2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： KRT8

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:

This gene is a member of the type II keratin family clustered on the long arm of chromosome 12. Type I and type II keratins heteropolymerize to form intermediate-sized filaments in the cytoplasm of epithelial cells. The product of this gene typically dimerizes with keratin 18 to form an intermediate filament in simple single-layered epithelial cells. This protein plays a role in maintaining cellular structural integrity and also functions in signal transduction and cellular differentiation. Mutations in this gene cause cryptogenic cirrhosis. Alternatively spliced transcript variants have been found for this gene.

Applications:	WB
Name of antibody:	KRT8
Immunogen:	Fusion protein of human KRT8
Full name:	keratin 8
Synonyms:	K8; KO; CK8; CK-8; CYK8; K2C8; CARD2
SwissProt:	P05787
WB Predicted band size:	54 kDa
WB Positive control:	Hela and A549 cell lysates
WB Recommended dilution:	10000-50000

