

## 小鼠抗 KRT18 单克隆抗体

中文名称： 小鼠抗 KRT18 单克隆抗体

英文名称： Anti-KRT18 mouse monoclonal antibody

别名： keratin 18; K18; CK-18; CYK18

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： KRT18

反应种属： Human, Monkey

标记物： Unconjugate

克隆类型： Mouse Monoclonal

### 技术规格

<b>Background:</b>	KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its filament partner keratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene.
<b>Applications:</b>	WB, IF
<b>Name of antibody:</b>	KRT18
<b>Immunogen:</b>	Fusion protein of human KRT18

<b>Full name:</b>	keratin 18
<b>Synonyms:</b>	K18; CK-18; CYK18
<b>SwissProt:</b>	P05783
<b>WB Predicted band size:</b>	48 kDa
<b>WB Positive control:</b>	COS7, 293T and Hela cell lysates
<b>WB Recommended dilution:</b>	1000-5000
<b>IF positive control:</b>	Hela cells
<b>IF Recommend dilution:</b>	500-800



