

## 小鼠抗 L3MBTL3 单克隆抗体

中文名称： 小鼠抗 L3MBTL3 单克隆抗体

英文名称： Anti-L3MBTL3 mouse monoclonal antibody

别名： L3MBTL3, histone methyl-lysine binding protein; MBT1; MBT-1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： L3MBTL3

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the malignant brain tumor (MBT) family of chromatin interacting transcriptional repressors. Members of this family function as methyl-lysine readers, which recognize methylated lysine residues on histone protein tails, and are associated with the repression of gene expression. The encoded protein may regulate hematopoiesis. Homozygous deletion of this gene has been observed in human patients with medulloblastoma.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	L3MBTL3
<b>Immunogen:</b>	Fusion protein of human L3MBTL3

<b>Full name:</b>	L3MBTL3, histone methyl-lysine binding protein
<b>Synonyms:</b>	MBT1; MBT-1
<b>SwissProt:</b>	Q96JM7
<b>WB Predicted band size:</b>	88 KD
<b>WB Positive control:</b>	Hela, 293T, A549, K562 cell lysates
<b>WB Recommended dilution:</b>	200-1000