

小鼠抗 LAMP1 单克隆抗体

中文名称：小鼠抗 LAMP1 单克隆抗体

英文名称：Anti-LAMP1 mouse monoclonal antibody

别名：lysosomal associated membrane protein 1; LAMPA; CD107a; LGP120

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：LAMP1

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.
Applications:	WB, FC
Name of antibody:	LAMP1
Immunogen:	Fusion protein of human LAMP1
Full name:	lysosomal associated membrane protein 1
Synonyms:	LAMPA; CD107a; LGP120
SwissProt:	P11279
WB Predicted band size:	45 KD



WB Positive control:	A431, COS7, 293T cell lysates
WB Recommended dilution:	200-1000
FC Recommend dilution:	50-100