

小鼠抗 LGR6 单克隆抗体

中文名称: 小鼠抗 LGR6 单克隆抗体

英文名称: Anti-LGR6 mouse monoclonal antibody

别名: GPCR; VTS20631

储存: 冷冻 (-20℃) 避光

抗原: LGR6

宿主: Mouse

反应种属: Human

标记物: Unconjugate

相关类别: 一抗

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes a member of the leucine-rich repeat-containing subgroup of the G protein-coupled 7-transmembrane protein superfamily. The encoded protein is a glycoprotein hormone receptor with a large N-terminal extracellular domain that contains leucine-rich repeats important for the formation of a horseshoe-shaped interaction motif for ligand binding. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2008].
Applications:	WB, IF
Name of antibody:	LGR6
Immunogen:	Fusion protein of human LGR6

Full name:	leucine-rich repeat containing G protein-coupled receptor 6 (LGR6), transcript variant 1
Synonyms:	GPCR; VTS20631
SwissProt:	Q9HBX8
WB Predicted band size:	104 kDa
WB Positive control:	Human uterus tissue lysate
WB Recommended dilution:	500-2000