

## 兔抗 TLE5 多克隆抗体

- 中文名称：兔抗 TLE5 多克隆抗体
- 英文名称：Anti-TLE5 rabbit polyclonal antibody
- 别名：TLE family member 5, transcriptional modulator; AE
- 相关类别：一抗
- 储存：冷冻 (-20℃)
- 抗原：TLE5
- 宿主：Rabbit
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is similar in sequence to the amino terminus of Drosophila enhancer of split groucho, a protein involved in neurogenesis during embryonic development. The encoded protein, which belongs to the groucho/TLE family of proteins, can function as a homooligomer or as a heterooligomer with other family members to dominantly repress the expression of other family member genes. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TLE5

<b>Immunogen:</b>	Fusion protein of human TLE5
<b>Full name:</b>	TLE family member 5, transcriptional modulator
<b>Synonyms:</b>	AES; GRG; ESP1; GRG5; AES-1; AES-2; Grg-5
<b>SwissProt:</b>	Q08117
<b>IHC positive control:</b>	Human liver cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	25-100

