

小鼠抗 LHX4 单克隆抗体

中文名称： 小鼠抗 LHX4 单克隆抗体

英文名称： Anti-LHX4 mouse monoclonal antibody

别名： LIM homeobox 4; CPHD4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： LHX4

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a member of a large protein family which contains the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein is a transcription factor involved in the control of differentiation and development of the pituitary gland. Mutations in this gene cause combined pituitary hormone deficiency 4.
Applications:	WB
Name of antibody:	LHX4
Immunogen:	Fusion protein of human LHX4
Full name:	LIM homeobox 4
Synonyms:	CPHD4

SwissProt:	Q969G2
WB Predicted band size:	43 KD
WB Positive control:	HepG2 cell lysate
WB Recommended dilution:	200-1000