

兔抗 TLE6 多克隆抗体

中文名称：兔抗 TLE6 多克隆抗体

英文名称： Anti-TLE6 rabbit polyclonal antibody

别名： transducin like enhancer of split 6; GRG6; PREMBL

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TLE6

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of the Groucho/ transducin-like Enhancer of split family of transcriptional co-repressors. The encoded protein is a component of the mammalian subcortical maternal complex, which is required for preimplantation development. In mouse, knock out of this gene results in cleavage-stage embryonic arrest resulting from defective cytoplasmic F-actin meshwork formation and asymmetric cell division. In human, an allelic variant in this gene is associated with preimplantation embryonic lethality. Alternative splicing results in multiple transcript variants.

Applications:	ELISA, IHC
Name of antibody:	TLE6
Immunogen:	Fusion protein of human TLE6
Full name:	transducin like enhancer of split 6
Synonyms:	GRG6; PREMBL
SwissProt:	Q9H808
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

