

小鼠抗 LIN28A 单克隆抗体

中文名称： 小鼠抗 LIN28A 单克隆抗体

英文名称： Anti-LIN28A mouse monoclonal antibody

别名： lin-28 homolog A; CSDD1; LIN28; LIN-28; ZCCHC1; lin-28A

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： LIN28A

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:

This gene encodes a LIN-28 family RNA-binding protein that acts as a posttranscriptional regulator of genes involved in developmental timing and self-renewal in embryonic stem cells. The encoded protein functions through direct interaction with target mRNAs and by disrupting the maturation of certain miRNAs involved in embryonic development. This protein prevents the terminal processing of the LET7 family of microRNAs which are major regulators of cellular growth and differentiation. Aberrant expression of this gene is associated with cancer progression in multiple tissues.

Applications:	WB
Name of antibody:	LIN28A
Immunogen:	Fusion protein of human LIN28A
Full name:	lin-28 homolog A
Synonyms:	CSDD1; LIN28; LIN-28; ZCCHC1; lin-28A
SwissProt:	Q9H9Z2
WB Predicted band size:	23 kDa
WB Positive control:	NTERA and F9 cell lysates
WB Recommended dilution:	1000-5000

